



Maybridge Chiral Resolution Screening and Purification Kits

Introduction

Diastereomeric crystallization is a commonly used effective process to obtain optically pure chiral compounds from their racemic mixtures. However, choosing the optimal conditions for the process; e.g., combination of resolving agents and solvents, is time-consuming, tedious and labor-intensive.

Maybridge Chiral Resolution Screening and Purification kits provide scientists with a quick and systematic approach to find the best separation conditions under which the target compound can be isolated with the highest yield and optical purity.

Key features and benefits

- Rapid Screening the kits include 384 different combinations of resolving agents and solvents, increasing the chances of finding the optimal separation conditions
- High Performance development time reduced to one day
- Efficient as little as 0.4mmol of racemate required
- **Ready to Use** resolving agents and solvents are pre-dispensed in 96-well plates
- Convenient the screening kits provide positive results identifiable by a quick visual or optical inspection, and the purification and recovery kit allows easy recovery and purification of the enantiomers



Types of Chiral Resolution Screening and Purification Kits

Product name	Description	Selection guide	Amount racemate required	Plate type	Product code
Maybridge Chiral	 4 x 96 plates containing 32 different acidic resolving agents 4 x 96 plates containing 12 different solvents 	Identifies optimal conditions for resolving base chiral compounds (e.g., amines)	12.0mmol	Plastic	CSA05004P
Resolution Screening Kit for base racemates			0.4mmol	Glass	CSA05004G
Maybridge Chiral Resolution Screening	• 2 x 96 plates containing 16 different strong	Identifies optimal conditions for resolving weak base chiral compounds (e.g., amines)	12.0mmol	Plastic	CSA05004STRG
Kit for weak base Racemates	acidic resolving agents • 2 x 96 plates containing 12 different solvents		0.4mmol	Glass	CSA05004GSTRG
Maybridge Chiral	 4 x 96 plates containing 32 different basic resolving agents 4 x 96 plates containing 12 different solvents 	Identifies optimal conditions for resolving acid chiral compounds (e.g., carboxylic acids)	12.0mmol	Plastic	CSB05004P
Resolution Screening Kit for acid racemates			0.4mmol	Glass	CSB05004G
Maybridge Chiral Purification and Recovery Kit for base racemates	1 x 200mL recovery solution for recovery of the resolved racemates 2 x 250mL bottles for filtration and purification 1 filter funnel 1 pack of pH paper	Used in conjunction with the screening kits, this kit provides the components necessary to recover and purify the chiral isomers from base racemates			CSA08004B
Maybridge Chiral Purification and Recovery Kit for acid racemates	1 x 200mL recovery solution for recovery of the resolved racemates 2 x 250mL bottles for filtration and purification 1 filter funnel 1 pack of pH paper	Used in conjunction with the screening kits, this kit provides the components necessary to recover and purify the chiral isomers from acid racemates			CSA08004A

How the kits work

Resolution of racemic chiral acids and amines can be achieved via the crystallization of their corresponding diastereomeric salt. The racemic compound, in this example a primary amine, is treated with an enantiomerically pure organic acid (in this case a carboxylic acid but sulfonic acids can be used as well).

Diastereomeric salt formation
$$\begin{array}{c} & & & & \\$$


Under the right conditions, one of the diastereomeric formed salts will precipitate as crystals whereas the other remains in solution. It is then possible to separate the salts by filtering the precipitate and recrystallizing from a suitable solvent. Further recrystallizations are often required to achieve the desired purity of the product.

The amine can then be isolated using an alkaline aqueous work-up to give the free base.

How to use the kits

- Dissolve the racemic starting material (12.0mmol required for the plastic vials and 0.4mmol for the glass vials) in a transfer solvent and distribute among the individual vials containing the resolving agents.
- Remove the transfer solvent via evaporation.
- Transfer the crystallization solvents to vials containing the racemic starting material and resolving agents.
- Heat the entire plate to ~80C for up to 5 minutes until a homogenous solution is achieved. (The plate is chemically inert and able to withstand extreme temperatures.)
- Cool the solutions until crystals begin to appear (either room temperature or in an ice/water bath).
- Analyze the crystals to identify which combination of resolving agent and solvent is optimal.
- At any point, recover the racemic material and extract and purify individual isomers using the components provided in the Purification and Recovery Kit.

Illustrated workflow





Austria

Belgium

For more information please visit our web site www.maybridge.com or contact your local distributor.

Denmark

Fisher Scientific (Austria) GmbH Rudolf von Alt-Platz 1 A-1030 Wien Toll-Free Number: 0800 20 88 40 Fax: 0800 20 66 90 info.austria@thermofisher.com www.at.fishersci.com	Fisher Scientific Clintonpark – Keppekouter Ninovesteenweg 198 B-9320 Erembodegem Tel: +32 (0) 56 260 260 Fax: +32 (0) 56 260 270 be.fisher@thermofisher.com www.be.fishersci.com	Fisher Scientific, spol. s r.o. Kosmonautu 324 CZ-530 09 Pardubice Tel: +420 466 798 230 Fax: +420 466 435 008 info.cz@thermofisher.com www.thermofisher.cz	Fisher Scientific Industrivej 3 DK-3550 Slangerup Tel: +45 70 27 99 20 Fax: +45 70 27 99 29 kundeservice@thermofisher.com www.fishersci.dk	Fisher Scientific Oy Ratastie 2 FI-01620 Vantaa Tel: +358 802 76 280 Fax: +358 802 76 235 fisher.fi@thermofisher.com www.fishersci.fi	Fisher Scientific FAS Parc d'Innovation BP 50111 67403 Illikirch Cedex Tel: +33 (0) 3 88 67 53 20 Fax: +33 (0) 3 88 67 11 68 fr.commande@thermofisher.com www.fr.fishersci.com
Germany	Ireland	Italy	The Netherlands	Norway	Portugal
Fisher Scientific GmbH Im Heiligen Feld 17 D-58239 Schwerte Toll-Free Number: 0800 3 47 43 70 Fax: 0800 3 47 43 71 info.germany@thermofisher.com www.de.fishersci.com	Fisher Scientific Ireland Suite 3 Plaza 212 Blanchardstown Corporate Park 2 Ballycoolin Dublin 15 Tel: +353 01 885 5854 Fax: +353 01 899 1855 fsie.sales@thermofisher.com www.ie.fishersci.com	Fisher Scientific Tel: +39 02 953 28 258 Fax: +39 02 953 27 374 it.fisher@thermofisher.com www.it.fishersci.com	Fisher Scientific Postbus 4 Scheepsbouwersweg 1b 1120 AA Landsmeer Tel: +31 (0) 20 487 70 70 Fax: +31 (0) 20 487 70 70 nl.info@thermofisher.com www.fishersci.nl	Fisher Scientific Frysjaveien 33E 0884 Oslo Tel: +47 22 95 59 59 Fax: +47 22 95 59 40 fisher.no@thermofisher.com www.fishersci.no	Fisher Scientific Rua Pedro Álvares Cabral, nº24, 3'D Edifício Euro - Infantado 2670-391 Loures Tel: +351 21 425 33 50/4 Fax: +351 21 425 33 51 pt.fisher@thermofisher.com www.pt.fishersci.com
Slovenia	Spain	Sweden	Switzerland	United Kingdom	
Fisher Scientific d.o.o. Parmova 53 SI-1000 Ljubljana Tel: +386 597 31 494 Fax: +386 597 31 495 si.info@thermofisher.com www.si.fishersci.com	Fisher Scientific C/ Luis I, 9 28031 Madrid Tel: +91 380 67 10 Fax: +91 380 85 02 es.fisher@thermofisher.com www.es.fishersci.com	Fisher Scientific Box 9193 400 94 Göteborg Tel: +46 31 352 32 00 Fax: +46 31 352 32 50 fisher.se@thermofisher.com www.fishersci.se	Fisher Scientific AG Wilstrasse 57 CH-5610 Wohlen Tel: +41 (0) 56 618 41 11 Fax: +41 (0) 56 618 41 41 info.ch@thermofisher.com www.ch.fishersci.com	Fisher Scientific UK Ltd. Bishop Meadow Road Loughborough Leicestershire LE11 5RG Tel: +44 (0)1509 555500 Fax: +44 (0)1509 555111 fsuk.sales@thermofisher.com www.fisher.co.uk	



Czech Republic



Finland

France

