

# New 25mL AcroSeal packaging

The recent introduction of a 25mL pack size to our industry-leading AcroSeal™ packaging range meets your needs for applications where small research-scale quantities of extra dry solvents are required.

Designed for your convenience, each bottle is packed using our new septum to ensure optimum performance in your moisture-sensitive work. The new 25mL package size expands our AcroSeal portfolio to better meet your chemistry needs.

## Highlights

- **Convenience:**  
Ideal for use when requiring small quantities
- **Easy to use:**  
Solvent is easily extracted with a standard needle
- **Quality:**  
Water content guaranteed at less than 50 parts per million



## Discover the products available in our new 25mL AcroSeal packaging

Product code	Product name	Pack size	CAS no.
448391000	Acetonitrile, 99.9%, Extra Dry	4 x 25mL	75-05-8
448371000	Dichloromethane, 99.9%, Extra Dry, stabilized	4 x 25mL	75-09-2
448421000	Diethyl ether, 99.5%, Extra Dry, stabilized	4 x 25mL	60-29-7
448381000	N,N-Dimethylformamide, 99.8%, Extra Dry	4 x 25mL	68-12-2
450171000	1,4-Dioxane, 99.8%, Extra Dry, stabilized	4 x 25mL	123-91-1
448411000	Methanol, 99.9%, Extra Dry	4 x 25mL	67-56-1
448361000	Tetrahydrofuran, 99.85%, Extra Dry, stabilized	4 x 25mL	109-99-9
448401000	Toluene, 99.85%, Extra Dry	4 x 25mL	108-88-3

## How to open the 25mL AcroSeal septum

Peel back the metal cap



Pierce the seal with the needle to withdraw the required amount of solvent



To get the best results from the septum, we recommend the use of 18 to 21 gauge needles.

### Related resources




#### AcroSeal product range overview

This brochure lists 600+ AcroSeal products, including solvents, reagents in solution, organics, organometallics and deuterated solvents.

Download a PDF via [www.acros.com](http://www.acros.com)

or request a printed copy via [perbio.euomarketing@thermofisher.com](mailto:perbio.euomarketing@thermofisher.com)



#### AcroSeal video

This video highlights how you can use our unique AcroSeal packaging to keep your air- and moisture-sensitive solvents and reagents drier for longer.

View the video at [www.acros.com](http://www.acros.com)

**ACRÖS**  
ORGANICS

©2013 Thermo Fisher Scientific Inc. All rights reserved.

These brands and trademarks are part of Thermo Fisher Scientific Inc. and its subsidiaries.

**Austria**  
info.austria@thermofisher.com  
www.at.fishersci.com  
Tel: 0800 20 88 40  
Fax: 0800 20 66 90

**Belgium**  
be.fisher@thermofisher.com  
www.be.fishersci.com  
Tel: 056 260 260  
Fax: 056 260 270

**Czech Republic**  
info.cz@thermofisher.com  
www.thermofisher.cz  
Tel: 466 798 230  
Fax: 466 435 008

**Denmark**  
kundeservice@thermofisher.com  
www.fishersci.dk  
Tel: 70 27 99 20  
Fax: 70 27 99 29

**Finland**  
fisher.fi@thermofisher.com  
www.fishersci.fi  
Tel: 09-802 76 280  
Fax: 09-802 76 235

**France**  
fr.commande@thermofisher.com  
www.fr.fishersci.com  
Tel: 03 88 67 53 20  
Fax: 03 88 67 11 68

**Germany**  
info.germany@thermofisher.com  
www.de.fishersci.com  
Tel: 0180 5 258221  
Fax: 0180 5 258223

**Ireland**  
fsie.sales@thermofisher.com  
www.ie.fishersci.com  
Tel: 01 885 5854  
Fax: 01 899 1855

**Italy**  
it.fisher@thermofisher.com  
www.it.fishersci.com  
Tel: 02 953 28 258  
Fax: 02 953 27 374

**The Netherlands**  
nl.info@thermofisher.com  
www.fishersci.nl  
Tel: 020 487 70 00  
Fax: 020 487 70 70

**Norway**  
fisher.no@thermofisher.com  
www.fishersci.no  
Tel: 22 95 59 59  
Fax: 22 95 59 40

**Portugal**  
pt.fisher@thermofisher.com  
www.pt.fishersci.com  
Tel: 21 425 33 50  
Fax: 21 425 33 51

**Slovenia**  
si.info@thermofisher.com  
www.si.fishersci.com  
Tel: + 49 180 5 258221  
Fax: + 49 180 5 258223

**Spain**  
es.fisher@thermofisher.com  
www.es.fishersci.com  
Tel: 902 239 303  
Fax: 902 239 404

**Sweden**  
fisher.se@thermofisher.com  
www.fishersci.se  
Tel: 31 352 32 00  
Fax: 31 352 32 50

**Switzerland**  
info.ch@thermofisher.com  
www.ch.fishersci.com  
Tel: 056 618 41 11  
Fax: 056 618 41 41

**United Kingdom**  
fsuk.sales@thermofisher.com  
www.fisher.co.uk  
Tel: 01509 555500  
Fax: 01509 555111



Part of Thermo Fisher Scientific