

Centrifuges

Fisherbrand bench top

Centrifuge, horizontal, Clinispin Horizon 842VES, Fisherbrand®

Fisherbrand®



For use by physicians, clinics, hospitals, veterinary clinics, fertility clinics, educational laboratories, industrial laboratories, reference laboratories and more.

- Performance Plus horizontal rotor and 5mL buckets included as standard
- Rapid separations
- Full horizontal separation
- Precise variable controls
- Easy operation
- Maintenance free induction motor
- Quiet operation and low vibration
- 10 memory settings

The Clinispin Horizon 842VES is the most compact of the performance series of centrifuges. The 842VES offers control over almost every aspect of the units operation; variable acceleration/deceleration, time and speed. The high speed and g-force capability produces platelet poor/free plasma in just two minutes and greatly reduces processing times. The unique horizontal rotor design allows for samples to be inserted at a convenient 45° angle, they are then gently rotated to a full horizontal position. The horizontal separation produces a better separation with a tighter gel line, eliminating remixing and re-spins.

Technical Specification

Capacity, mL	6 x 15
Speed, rpm	6,500
RCF, x g [max.]	5,000
Timer range, min	1 to 99
Dimensions (w x d x h), mm	280 x 340 x 200
Mass, kg	10.5
Electrical supply	240V, 50/60Hz

Rotor	Tube size [max.]	Speed, max., rpm	RCF, x g [max.]
6 place Performance Plus horizontal rotor	5mL, 17mm x 75mm	6,500	5,000
6 place horizontal rotor	10mL, 17mm x 10mm	4,500	2,800
6 place fixed angle rotor	15mL, 17mm x 125mm	4,000	1,950

Catalogue No	Legacy No	Description
12375407	-	Clinispin Horizon 842VES centrifuge

Accessories

Catalogue No	Legacy No	Description
12385407	-	Performance Plus horizontal rotor, 6 place
12395407	-	Bucket, 15mL (orange), for maximum tube size 17mm x 125mm
12305417	-	Bucket, 10mL (purple), for maximum tube size 17mm x 100mm
12315417	-	Bucket, 5mL (blue), for maximum tube size 17mm x 75mm
12325417	-	Bucket, 5mL Performance Plus, for maximum tube size 17mm x 75mm