



Fisherbrand[®]
QUALITY. RELIABILITY. VALUE.

Glassware

Disposable, Reusable and Specialty



ThermoFisher
SCIENTIFIC

Fisherbrand Class 'A' Volumetric ware – Value added

- ▶ Individual Calibrated & Certified with distinct Identification number
- ▶ Conforming to international Specifications ISO / ASTM
- ▶ Made from Precision Bore Tubing
- ▶ Graduation in Amber
- ▶ Complete Traceability



ISO

ASTM / USP

Wide Range of Burettes & Pipettes

Burettes: 10 mL - 100 mL
Straight bore, Rota flow, Glass Key, PTFE Key

Pipettes: 1 mL - 100 mL
Complete Delivery, Mohr, One mark / Bulb



Burettes Class 'A'		
Capacity (mL)	Graduation (mL)	Tolerance (± mL)
10	0.05	0.02
25	0.10	0.05
50	0.10	0.05
100	0.20	0.10

Pipettes Class 'A'		
Capacity (mL)	Graduation (mL)	Tolerance (± mL)
1	0.01	0.006
2	0.02	0.010
5	0.05	0.030
10	0.10	0.050
25	0.20	0.100

Transfer Pipettes Class 'A'	
Capacity (mL)	Tolerance (+ mL)
1	0.007
2	0.010
5	0.015
10	0.020
20	0.030
25	0.030
50	0.50
100	0.100

Volumetric Flasks Class 'A'

Clear & Amber : 5 mL - 2000 mL

Capacity (mL)	Tolerance (± mL)
1	0.015
2	0.015
3	0.015
5	0.025
10	0.025
20	0.04
25	0.04
50	0.06
100	0.10
200	0.15
250	0.15
500	0.25
1000	0.40
2000	0.60



Measuring Cylinders Class 'A'

Hexagonal base, Round base with I / C Stopper
10 mL - 1000 mL

Capacity (mL)	Graduation (mL)	Tolerance (± mL)
10	0.2	0.10
25	0.5	0.25
50	1.0	0.50
100	1.0	0.50
250	2.0	1.00
500	5.0	2.50
1000	10.0	5.00



Fisherbrand Boilingware

An Exhaustive Range...

- Uniform wall Thickness
- Thermal Shock Proof
- Made from 3.3 Borosilicare Glass


Flasks

Round Bottom, Flat Bottom,
Conical with & without I/C Joints
Multiple Necks - Parallel & Angled



Condensers

Liebig, Allihn & Graham



Soxhlet

Soxhlet Extraction Apparatus
Distillation Apparatus - Ammonia



Products for Special Applications

- Desiccators Vacuum Top & Knob Top
- Flask for suger Estimation
- Haffkine Flasks 3000 mL - 4000 mL
- Iodine Flasks 250 mL - 500 mL
- Kjeldahl Flasks 500 mL - 800 mL
- Museum Jars Wide Range / Custom Sizes
- Seprating Funnels Pear shape PTEF /glass key, Straight / Rota Flow
- Filtration flasks



Products for General Applications


- Beakers
- Conical Flasks
- Cell Culture Tubes with Screw Cap, Teflon Wad, Flat Bottom, Round Bottom
5 mL - 60 mL
- Wide Range of Centrifuge Tubes
15 mL - 50 mL
- Test Tubes with Rim & without Rim
- Stoppered Test Tubes Plain & Graduated
5 mL - 100 mL
- Petri dishes autoclavable
90 x 15 mm, 100 x 15 mm & 100 x 20 mm



Container ware – Fisherbrand bottles

Mobile Phase Reservoirs

- Volume Range 25 mL - 20000 mL
- Clear & Amber
- Enhanced graduations
- Permanent white enamel marking spot
- Optically clear, drip-free
- replaceable polypropylene pour ring
- Linerless, autoclavable (to 140°C) polypropylene cap



Fisherbrand Bottles

Wide Range of Bottles to suit varied Application

- Narrow mouth Reagent Bottles, Plain & Amber
60 mL – 2000 mL
- Weighing Bottles, Relative Density Bottles
- BOD Bottles – 125 mL & 300 mL
- Screw Cap with serial number, Wide mouth, Reagent Bottles, Plain & Amber
100 mL – 20000 mL
- Tooled Neck Bottles
250 mL -20000 mL
- Aspirator with Stop Clock
5000 mL – 20000 mL

Filtration Assembly

- Excellent Optical Clarity provides
- Provides more effective Filtration Area
- Autoclavable
- Heavy Base
- Base Study Design

Price on request

Solutions...Unlimited...
For more details & exciting offers contact



© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc., and its subsidiaries.

Thermo Fisher Scientific India Pvt. Ltd.

403-404, Delphi 'B' Wing, Hiranandani Business Park, Powai, Mumbai - 400 076.

Phone No.: 022-6680 3000 • Fax: 022-6680 3001/2

Indi Toll Free No.: 1-800-222-230

Email : lsi_marketing@thermofisher.com

EMFBG62014FY

ThermoFisher
SCIENTIFIC