

Fisherbrand* Standard Mini-Centrifuge

The Fisherbrand mini-centrifuge is ideal for microfilter cell separations and HPLC samples. Operation is now faster and easier and the ergonomically designed lid will fit in the palm of the hand for easy repeat opening and use.

Includes:

- Two easy change, tool-free rotors: Six-place standard rotor for 1.5/2mL tubes, Eight-place strip rotor for two 0.2mL tube strips or 16 single tubes
- Six 0.2mL and six 0.5mL adapters for customized use
- Storage Case — Keeps rotors and adapters contained and dust-free
- Tube Rack — Convenient for handling samples before and after quick spins



Description
<ul style="list-style-type: none"> ▪ Features exclusive quick-release rotor for fast and easy changes ▪ Internal safety switch ensures unit will not operate without lid in place ▪ Removable cord for compact storage ▪ Maximum speed: 6000rpm ▪ Max RCF: 2000G

Cat. No.	Maximum Speed	Description
13406188	6000rpm	Standard Mini-Centrifuge Mini-centrifuge includes 6 place rotor for 1.5mL/2.0mL tubes & 8 place rotor for 2 tube strip UV resistant plastic 220-230V 50/60Hz 6,000rpm maximum speed 153mm x 128mm x 104mm Fisherbrand