



Laboratory Carts & Benches



*Protecting your
laboratory environment*

LABCONCO

Labconco Laboratory Carts and Benches fit in every lab.

These versatile carts have endless uses from transporting equipment and supplies to providing auxiliary bench space.


Labconco Carts and Benches feature:

- Ready-to-roll, full factory assembly.
- Toe-locking casters.
- Ergonomic, molded plastic handle grips.
- Durable and chemical-resistant, glacier white, epoxy-coated tubular steel frames.
- ID tag to ensure your cart makes a round trip back to you.
- Lifetime guarantee – if your Labconco Cart or Bench does not live up to your expectations, we will replace it FREE.

Catalog #	Description	Shipping Weight
8010000	Flexi-Bin Cart with 4" bin	55 lbs. (25 kg)
8010001	Flexi-Bin Cart with 10" bin	55 lbs. (25 kg)
8068000	Bin, 4" deep, with dividers	5 lbs. (2.3 kg)
8068100	Bin, 10" deep, with dividers	10 lbs. (4.5 kg)



#8010000



#8010001

Flexi-Bin™ Carts

- For safe transport of glassware and chemical containers.
- Removable, chemical-resistant polyethylene bin, 22.25 x 17.0" (56.5 x 43.2 cm), is included. Inside dimensions: 20.25 x 15.25" (51.4 x 38.7 cm). Bins are offered in two overall depths of 4" (10 cm) and 10" (25 cm) with inside depths of 3.8" (9.7 cm) and 9.8" (24.9 cm), and are interchangeable.
- Include three long (20.25") and four short (15.25") polyethylene dividers for partitioning the bin to securely hold items in upright positions.
- Bottom open shelf holds additional supplies.
- 4" diameter (10 cm) casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.9 x 19.0 x 36.4" high (88.6 x 48.3 x 92.4 cm)

See dimensional data on page 8.

Durable

The 18-gauge welded steel frames of 1" square steel tubing are painted with tough, corrosion-resistant, dry powder epoxy.

Benches have 2" square steel tubing.

Stockroom Cart

- Catalog #8021000.
- Step ladder and cart in one.
- Two removable, watertight, epoxy-coated steel pans with gray vinyl mats, and one vinyl-coated wire basket.
- Sturdy, two-step ladder with non-skid surface conforms to OSHA and ANSI design guidelines.
- Spring-loaded, 9.5" (24.1 cm) deep steps with skid-resistant, rubber feet prevent cart from moving while user stands on the ladder.
- Steps withstand loads up to 1000 lbs. (453.6 kg).
- 4" diameter (10 cm) casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: Cart is 19.0 x 39.5 x 43.0" high (48.3 x 100.3 x 109.2 cm). Top pan is 6.5 x 16.75 x 2.0" (16.5 x 42.5 x 5.1 cm). Bottom pan is 17.5 x 16.75 x 2.0" (44.5 x 42.5 x 5.1 cm). Basket is 16.0 x 18.0 x 4.0" (40.6 x 45.7 x 10.2 cm).
- Shipping weight 70 lbs. (32 kg).



#8021000



#8007000

Utility Cart

- Catalog #8007000.
- Upper open shelf transports and supports small electrical equipment. May be used as a compact work station for laboratory equipment.
- Lower 2.0" (5.1 cm) deep epoxy-coated steel, removable watertight pan transports wet or soiled glassware and other labware.
- Power strip with 115 volts, 60 Hz, 12 amps capacity, has 6 electrical outlets, cord bracket and a 12-foot, 3-wire cord with plug.
- Two corrosion-resistant, one-piece, vacuum-formed ABS plastic drawers, 15.75 x 14.9 x 4.2" high (40.0 x 37.8 x 10.7 cm), hold small labware and supplies.
- 4" diameter (10 cm) casters, two with toe locks.
- Supports loads up to 400 lbs. (181 kg).
- Lifetime guarantee.
- Dimensions: 34.9 x 19.0 x 36.4" high (88.6 x 48.3 x 92.4 cm)
- Shipping weight 70 lbs. (32 kg)

See dimensional data on page 8.

Portable Table

- Catalog #8025000.
- Labconco's most popular cart has two open shelves that provide easy transfer of equipment from countertop to cart.
- 4" diameter (10 cm) casters, two with toe locks.
- Supports loads up to 400 lbs. (181 kg).
- Lifetime guarantee.
- Dimensions: 34.9 x 19.0 x 36.4" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 50 lbs. (23 kg).



#8025000



#8020000

Chemical Cart

- Catalog #8020000.
- Similar to the Portable Table except with 1.2" (3.0 cm) high rigid PVC extrusions on two sides to prevent containers from sliding off. The lower shelf also has rigid plastic guard rails on both ends.
- 4" diameter (10 cm) casters, two with toe locks.
- Lifetime guarantee.
- Supports loads up to 400 lbs. (181 kg).
- Dimensions: 34.9 x 19.0 x 36.4" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 50 lbs. (23 kg).

Stainless Steel Cart

- Catalog #8019000.
- Holds a variety of laboratory supplies, equipment and miscellaneous containers.
- Two removable polished Type 304 stainless steel pans are 2" (51 mm) deep on three sides with one open side for easy equipment transfer from countertop to cart.
- 4" diameter (10 cm) casters, two with toe locks.
- Supports loads up to 400 lbs. (181 kg).
- Lifetime guarantee.
- Dimensions: 34.9 x 19.13 x 36.4" high (88.6 x 48.6 x 92.4 cm).
- Shipping weight 50 lbs. (23 kg).



#8019000

Quiet

Four-inch diameter casters absorb vibration over rough, uneven floors. Two toe-locks keep carts stationary when needed. Benches have 5" diameter, non-marking toe-locking casters with bearings for quiet operation.



Glassware Carts

- Carry glassware, plasticware, even small instruments, in vinyl-coated wire baskets.
- Choose 2-basket or 4-basket style.
- Complete with removable plastic drip pan to catch spills.
- 4" diameter (10 cm) casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 34.9 x 19.0 x 36.4" high (88.6 x 48.3 x 92.4cm).

Catalog #	Description	Shipping Weight
8032500	2-Basket Glassware Cart. Each basket measures 32.0 x 18.0 x 4.0" (81.3 x 45.7 x 10.2 cm)	50 lbs. (23 kg)
8045000	4-Basket Glassware Cart. Each basket measures 16.0 x 18.0 x 4.0" (40.6 x 45.7 x 10.2 cm)	50 lbs. (23 kg)
8040100	Large Basket, 32.0 x 18.0 x 4.0" (81.3 x 45.7 x 10.2 cm)	5 lbs. (2.3 kg)
8040200	Small Basket, 16.0 x 18.0 x 4.0" (40.6 x 45.7 x 10.2 cm)	2 lbs. (0.9 kg)

Instrument Desk

- Catalog #8052500.
- Serves as a desk, instrument table or laboratory cart.
- Support bar is centered to permit comfortable seating from either side.
- Work surface is phenolic board with a melamine high-pressure laminate surface that resists stains and corrosion.
- 4" diameter (10 cm) casters, two with toe locks.
- Lifetime guarantee.
- Dimensions: 42.13 x 19.00 x 31.63" high. (107 x 48.3 x 80.3 cm).
- Shipping weight 40 lbs. (18 kg).



#8052500



#8000300

Carboy Caddy

- Catalog #8000300.
- Safe way to move large carboys or cartons.
- Helps eliminate heavy lifting and back strain.
- 3" diameter (8 cm) casters provide low profile, elevating shelf only 4.25" (10.8 cm) above the floor.
- Shelf is made of phenolic board with a melamine high-pressure laminate surface.
- Handle measures 38.5" high and 19.0" wide (97.8 x 48.3 cm).
- Shelf dimensions: 20.0 x 23.6" (50.8 x 60 cm).
- Supports loads up to 360 lbs. (163.3 kg).
- Lifetime guarantee.
- Shipping weight 35 lbs. (16 kg).

Pan Cart

- Catalog #8047500.
- Two 2" deep (5 cm) epoxy-coated steel, removable watertight pans transport wet or soiled glassware and other labware.
- Removable gray rubber mats line pans, prevent loads from shifting.
- 4" diameter (10 cm) casters, two with toe locks.
- Supports loads up to 400 lbs. (181 kg).
- Lifetime guarantee.
- Dimensions: 34.88 x 19.00 x 36.38" high (88.6 x 48.3 x 92.4 cm).
- Shipping weight 50 lbs. (23 kg).



#8047500

Strong

Labconco Carts support heavy loads.
Carts hold up to 400 pounds; Benches up to 540.

Mobile Bench

- Catalog #8060000.
- Heavy duty auxiliary bench may be transported from lab to lab.
- Durable work surface is fiber board with high pressure melamine laminate.
- 5" diameter (13 cm), non-marking polyurethane casters with toe locks.
- Supports loads up to 540 pounds (245 kg).
- Lifetime guarantee.
- Dimensions: 48.0 x 24.0 x 37.2" high. (121.9 x 61 x 94.5 cm). Work surface is 1.1" (29 mm) thick.
- Shipping weight 105 lbs. (48 kg).



#8060000

Variable Height Bench

- Catalog #8075000.
- Heavy duty auxiliary bench has adjustable height.
- Durable work surface is fiber board with high pressure melamine laminate.
- 5" diameter (13 cm), non-marking polyurethane casters with toe locks.
- Four adjustments from 29.6 to 37.4" (75.1 to 94.9 cm) using steel pins with nuts inserted in vertical frame members.
- Supports loads up to 540 pounds (245 kg).
- Lifetime guarantee.
- Work surface dimensions: 38.0 x 28.0 x 1.1" thick (96.5 x 71.1 cm x 29 mm).
- Shipping weight 105 lbs. (48 kg).


See dimensional data on page 8.




#8075000

D I M E N S I O N A L D A T A


Flexi-Bin™ Carts - #8010000 & #8010001




Glassware Carts - #8045000 & #8032500




Stockroom Cart - #8021000




Instrument Desk - #8052500




Utility Cart - #8007000




Carboy Caddy - #8000300




Portable Table - #8025000




Pan Cart - #8047500




Chemical Cart - #8020000




Mobile Bench - #8060000



Stainless Steel Cart - #8019000



Variable Height Bench - #8075000




Labconco offers additional options to add mobility to your laboratory.

- **Mobile Equipment Tables** come ready to roll and support portable ventilation enclosures, instruments or other bulky apparatus up to 800 pounds.
- **Hydraulic Lift Base Stands** may be combined with **Flat Solid Epoxy Work Surfaces** and **Caster Kits** to achieve portable auxiliary bench space with the added benefit of infinite height adjustment as needed.
- The **Logic® Cart**, originally designed as an accessory for Purifier Logic Biosafety Cabinets, has endless uses in the laboratory. Its durable polyethylene and ABS construction holds up to many laboratory environments.
- The **ScrubberMate™ Glassware & Rack Cart** is an accessory for the **SteamScrubber®** and **FlaskScrubber® Glassware Washers**, but could also be used as a portable carrier for beakers and flasks.

Mobile Equipment Tables

- May be used to transport heavy equipment or to provide auxiliary bench space.
- Base is constructed of durable 1.75" tubular steel with a chemical-resistant epoxy coating.
- Attractive glacier white exterior complements laboratory casework and all Labconco products.
- Telescoping legs adjust to eight height positions in 1" increments from 28.75" to 35.75" (73.0 to 90.8 cm).
- 5" diameter non-marking polyurethane casters with toe locks have bearings for quiet operation.
- Supports loads up to 800 pounds.
- Work surface is seamless, chemical-resistant flat solid epoxy, 1.25" thick and 30" deep.
- Available in 4, 5 and 6' widths.
- Fully assembled. Ready to roll.




Mobile Equipment Tables			
Nominal Width	Catalog #	Overall Width x Depth	Shipping Weight
4 feet	3819100	48.0 x 30.0" (121.9 x 76.2 cm)	250 lbs. (113.4 kg)
5 feet	3819101	60.0 x 30.0" (152.4 x 76.2 cm)	300 lbs. (136.1 kg)
6 feet	3819102	72.0 x 30.0" (182.9 x 76.2 cm)	350 lbs. (158.8 kg)



Hydraulic Lift Base Stands

- When coupled with a work surface, the stands may be used as an adjustable height table. Work Surfaces are sold separately.
- Durable 1.75" (4.4 cm) tubular, epoxy-coated steel frame with rear crossbar supports. Sub-frame is 2.0 x 3.0" (5.1 x 7.6 cm).
- Working height is infinitely adjustable from 26.75 to 34.75" high (67.9 x 88.3 cm).
- Four powder-coated, die-cast feet with levelers.
- Supports loads up to 800 pounds (363 kg).
- Optional accessory 5" diameter (13 cm), toe-locking Caster Kit is available for field installation.
- ADA*-compliant.

3784000 Caster Kit. Includes four 5" diameter (13 cm), non-marking polyurethane casters with toe locks. Casters have bearings for quiet operation. Casters replace the fixed feet that are included with the stand. **Installation is required.** Shipping weight 10 lbs. (4.5 kg)



#3780400 shown with #4882805 and #3784000.



* Americans with Disabilities Act
** International electrical configuration

Flat Solid Epoxy Work Surfaces

Nominal Width	Catalog #	Overall Width x Depth x Thickness	Shipping Weight
3 feet	4882805	36.0 x 30.0 x 1.25" (91.4 x 76.2 x 3.2 cm)	100 lbs. (45 kg)
4 feet	4882803	48.0 x 30.0 x 1.25" (121.9 x 76.2 x 3.2 cm)	135 lbs. (61 kg)
6 feet	4882804	72.0 x 30.0 x 1.25" (182.9 x 76.2 x 3.2 cm)	200 lbs. (91 kg)

Manual Hydraulic Lift Base Stands


Nominal Width	Catalog #	Overall Width x Depth	Shipping Weight
3 feet	3780400	36.0 x 29.0" (91.4 x 73.7 cm)	150 lbs. (68 kg)
4 feet	3780401	48.0 x 29.0" (121.9 x 73.7 cm)	180 lbs. (82 kg)
6 feet	3780402	72.0 x 29.0" (182.9 x 73.7 cm)	220 lbs. (100 kg)

Electric Hydraulic Lift Base Stands


Nominal Width	Catalog #	Overall Width x Depth	Electrical Requirements	Shipping Weight
3 feet	3780300	36.0 x 29.0" (91.4 x 73.7 cm)	115 volts, 60 Hz, 2.4 amps	150 lbs. (68 kg)
3 feet	3780303**	36.0 x 29.0" (91.4 x 73.7 cm)	230 volts, 50 Hz, 0.8 amp	150 lbs. (68 kg)
4 feet	3780301	48.0 x 29.0" (121.9 x 73.7 cm)	115 volts, 60 Hz, 2.4 amps	180 lbs. (82 kg)
4 feet	3780304**	48.0 x 29.0" (121.9 x 73.7 cm)	230 volts, 50 Hz, 0.8 amp	180 lbs. (82 kg)
6 feet	3780302	72.0 x 29.0" (182.9 x 73.7 cm)	115 volts, 60 Hz, 2.4 amps	220 lbs. (100 kg)
6 feet	3780305**	72.0 x 29.0" (182.9 x 73.7 cm)	230 volts, 50 Hz, 0.8 amp	220 lbs. (100 kg)

D I M E N S I O N A L D A T A

Logic® Cart - #8022000



ScrubberMate™ Glassware & Rack Cart - #8027000





#8022000

Logic[®] Cart

- Catalog #8022000.
- Handy cart stores lab supplies and may be moved easily to within reach of the user working at a bench, fume hood or biological safety cabinet.
- Three right side compartments, 3.5 x 17.0 x 4.5" high (9 x 43 x 51 cm), accommodate pipettes and other supplies.
- Top left side compartment with front panel cutouts accommodates boxes of wipes and gloves. The hinged panel may be pivoted down for loading.
- Drawer, 14.5 x 15.5 x 4" high (37 x 39 x 10 cm), holds a variety of miscellaneous small instruments and supplies.
- Open storage area in bottom, 15.8 x 16.0 x 12.0" high (40 x 41 x 30 cm), holds larger equipment and supplies.
- Work surface and bottom are constructed of corrosion-resistant white high density polyethylene that is easy to clean. Remaining components are constructed of corrosion-resistant white ABS.
- Compact dimensions allow easy fit under Purifier Logic Biosafety Cabinets mounted on Labconco Base Stands as well as other biological safety cabinets on open bases.
- Work surface is 25.0 x 19.0" deep (64 x 48 cm). Cabinet is 21.5 x 17.4" (55 x 44 cm). Overall height is 27.0" (68.6 cm).
- Work surface has handle cutouts in front and right hand side for easy access.
- 2" diameter (5 cm) casters provide mobility.
- Fully assembled.
- Shipping weight 50 lbs. (23 kg)



Logic Carts are shown stored beneath two Purifier Logic Biosafety Cabinets on Hydraulic Lift Base Stands.

See dimensional data on page 10.



#8027000

ScrubberMate[™] Glassware & Rack Cart

- Catalog #8027000.
- Provides convenient storage and transport of extra Labconco washer racks, inserts or miscellaneous glassware.
- 48 adjustable and removable pins, 0.5" wide x 4.75" high (1.3 x 12.1 cm), may be arranged in the drawers to hold inverted flasks, beakers and other glassware. Alternatively, accommodates up to two SteamScrubber[™] or FlaskScrubber[™] Glassware Washer racks or four wire inserts.
- Rails are provided to allow a lower washer rack to slide directly onto the cart from the washer.
- The lower drawer pulls out for simple loading of inserts or individual glassware.
- Durable and chemical-resistant, high density polyethylene construction.
- Includes 2.5" diameter casters with toe locks.
- Dimensions: 23.3 x 22.6 x 31.6" high (59.2 x 57.4 x 80.3 cm).
- Fully assembled.
- Shipping weight 70 lbs. (32 kg).



Sixteen combs, each with three pins, may be placed in sixteen positions in the upper drawer and sixteen positions in the lower drawer.



Upper and lower drawers hold two inserts each and pull out for easy access. An upper rack may be stored on the top; a lower rack below.

Visit www.labconco.com or call Labconco at 1-800-821-5525 for additional information about these quality products for your laboratory.



Fume Hoods & Ductless Enclosures



HEPA-Filtered Safety Cabinets, Enclosures & Clean Benches



Balance & Bulk Powder Enclosures



Nanotechnology Enclosures



Laboratory Animal Research Stations



Forensic Enclosures



Glove Boxes



Freeze Dry Systems



Centrifugal Vacuum Concentrators & Cold Traps



Multiple Sample Evaporation Systems



Glassware Washers



Water Purification Systems



Agricultural Chemistry Products



Blood Drawing Chairs



Vacuum Desiccator



Labconco Corporation • 8811 Prospect Avenue • Kansas City, Missouri 64132-2696
800-821-5525 • 816-333-8811 • Fax: 816-363-0130 • labconco@labconco.com • www.labconco.com

© 2010 by Labconco Corporation

Printed in the U.S.A. Product design subject to change without notice.

2-58-9/10-4M-Mainline-R1